



**beyond the technology**



**CSols**  
laboratory systems integration

## **RFID Innovative Sample Labelling, Scanning and Tracking Technologies**

**CSols RFTrackIT – Capabilities & Demo**



# RFID Team

- Team:

- Keith Huxford

CSols RTrackIT Brand Manager

- Frits van Calker & Andy Ho

RFIDDirect MD & Problem Solver


- Cedric Mochon

CSols Project Manager, Developer

# Agenda


- Introduction / CSols Capabilities / Recent Projects / New Projects
- An introduction to the World of RFID
  - Bottle scanning demo
  - Shoe shop video – finding items
- RFTtackIT – crate scanning, expected and not expected bottles
- What next?

# The First Application of RFID Technology



# Radio Frequency Identification (RFID)

Uses electromagnetic waves to transfer information to and from 'tags' or 'smart labels'...



... which can be attached to or embedded within a wide range of objects or solutions

# RFID Tag Types

- **Passive**

- Operational power scavenged from reader radiated power



- **Semi-passive**

- Operational power provided by battery







- **Active**

- Operational power provided by battery - transmitter built into tag




# RFID – How it works



- Reader broadcasts signal through antenna to transponder
- Transponder is charged with enough energy to send back identifying response
- Reader sends info/data to computer system for logging and processing

# RFID Scanners





# RFID read / write label printers

## RFID Printer Technology

- High Speed Printers
  - Can print on paper RFID labels
  - Sample information, 2D barcodes etc.
- CSols software 'pairs' sample id with the RFID tags UID and stores the data in the LabCache database



# Applications

# CSols RFID Applications



- Crate & Bottle Scanning
- Chain of Custody – Courier transfers, Night Stores, inter van, lab delivery
  - Temperature monitoring of vans / samples
- People, Column, Instrument Tracking
  - Asset / High Value item tracking
  - What was last run on column and how?
- Gold Ore analysis / traceability
- Solution tagging and tracking
- MS Azure MRC UK Brain Bank specimen storage in ultra low temperature freezers and at room temperature
  - Retrospectively labelling frozen samples

# CSols Gold Lab Testing


- ✦ Track all sacks, bags, sub-samples



# Flexible RFID/Barcode Sampling Templates




# Biopharma / MRC UK Brain Bank Network



**Flexible Sampling Templates**  
Also control “Tri-Safe” labelling



Cryogenic RFID label



Room Temperature  
QR Code RFID label

# Water company demo configuration

Site Browser

Search

- Demo site
  - Laboratory
  - Night Store
  - Van 1 (clean)
    - Consignment\_Ref123
  - Van 2 (waste)

Sample Name	Internal Reference	Status
Clean Bottle00007918	Demo Reference #123	Available
Clean Bottle00007919	Demo Reference #123	Available
Clean Bottle00007920	Demo Reference #123	Available
Clean Bottle00007921	Demo Reference #123	Available
Clean Bottle00007924	Demo Reference #123	Available
Clean Bottle00007925	Demo Reference #123	Available

Sample Browser

Search | Clear | New | Export | Import | Toggle Multi Select

- Demo Reference #123\_Consignment00007917
  - Demo Reference #123\_Clean Bottle00007918
  - Demo Reference #123\_Clean Bottle00007919
  - Demo Reference #123\_Clean Bottle00007920
  - Demo Reference #123\_Clean Bottle00007921
  - Demo Reference #123\_Clean Bottle00007924
  - Demo Reference #123\_Clean Bottle00007925

Sample Name	Internal Reference	Location
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123
Clean Bottle000079...	Demo Reference #123	\Van 1 (clean)\Consignment_Ref123

# Consignment Tracking

- Flag too many too few samples in a consignment
- Flag Items in wrong place / incorrect samples in van compartment (waste/clean)
- Prevent Clean in waste van and vice versa before it happens
- Flag ad hocs in a consignment
- Prevent accidental discards
- Notify of required discard(s)

Allocate discard place

Bottle RFID Tag	Sample ID	status
300EFE2F94D0FB4000005A02	10068	Wanted
E2001021400C01952360235E	10069	Found
202020204E384B3643494547	10070	Found
00006806000000000005555	10071	Found
E2001021440902721180A2B4	10072	Found
E20010214409023411909DC9	10073	Wanted
E2001021440902541180A26C		Unknown
300833B2DD0906C000000000		Unknown
E2001021440902731180A2B8		Unknown



# Flexible RFID/Barcode Sampling Templates

Site Survey : Wan 1 (clean)/Consignment\_Ref123

Sample ID	Int Ref	Found	Location
Expected Samples for (Wan 1 (clean)/Consignment_Ref123)			
Clean Bottle0007918	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007919	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007920	Demo Reference #123	False	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007921	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007924	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007925	Demo Reference #123	False	Wan 1 (clean)/Consignment_Ref123

Wan 1 (clean)/Consignment\_Ref123 | Sub Locations Not Included | Unidentified : 0 | Unexpected : 0 | Found : 4 of 6

Two missing samples

Site Survey : Wan 1 (clean)/Consignment\_Ref123

Sample ID	Int Ref	Found	Location
Expected Samples for (Wan 1 (clean)/Consignment_Ref123)			
Clean Bottle0007918	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007919	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007920	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007921	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007924	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Clean Bottle0007925	Demo Reference #123	True	Wan 1 (clean)/Consignment_Ref123
Unexpected samples for (Wan 1 (clean)/Consignment_Ref123)			
Clean Bottle0007922	Demo Reference #123	False	Wan 1 (clean)/Consignment_Ref123

Wan 1 (clean)/Consignment\_Ref123 | Sub Locations Not Included | Unidentified : 0 | Unexpected : 1 | Found : 6 of 6 | Sample 7918 3034000000000B0C0000061CB

One wrong sample

# CSols Applications / Ideas



- Integrated Remote Sampling Room, Van, Fridge, Freezer temperature monitoring
- Chemicals & Reagents Tracking
- Instrument / Personnel / Asset Inventory Tracking
- Column Identification / Instrument Information / PAT
- Gas Cylinder Identification & Tracking
- Linking items together – Instrument Modules, Pipes / Connectors, People / Tasks, Fuel Delivery Systems
- Automated PrepLab login and weighing / instrument queue building and setup

# Temperature Monitored Consignments



## Fully Integrated Vertical Audit Reports



CSols Remote Sampler



# Temperature Monitored Consignments, Boxes, Bags etc.



# Novel Techniques – ‘Solution tagging’

## Challenges

- Identifying ore / rock cyanide slurry sample id



# Tagging Columns for better information

- Also tag instruments / assets, columns and components to ensure modular instrument is validated
- Feedback support, management and maintenance data in real time
- Feedback real time usage / service data across company



# NFC / HF Tags with stored data

73% 1:40 PM

## Instrument ID

LC-MS001A  
Model E1M18  
23 Jan 2013  
A343-24-904  
36kgs  
SN 0171698476  
Organics Lab  
Dave Smith  
01827 456322

82% 12:22 PM

## Personnel Authorization

Keith Huxford  
3636756  
keith.huxford@csols.com  
www.csols.com  
01928 246233  
Sales & Marketing



78% 12:55 PM

## Column ID

KH179415  
367  
M1344-34A  
ODS HYPERSIL 5  
100X2.1  
B4678A  
pH 12, 350 bar, 70c



NFC Scanning Area

# Another Project Idea

## Automated PrepLab

- Automatic user authentication
- Automated sample identification & weighing
- QC samples added
- Automated user worklist build & automated instrument queue build





# Key Benefits of RFID

- Carry unique identifiers that positively identify items
- Saves time, speeds up tasks, improves traceability
- Does not require direct line of sight to be read
- Multiple tags can be read at one time
- Data can be stored in tags
- Can easily track movement of multiple items between locations
- Can automate actions when tag is read (cloud links, SMS, alarms, email, link to results/invoice etc.)



beyond the technology



**CSols**  
laboratory systems integration

## RFID Innovative Sample Labelling, Scanning and Tracking Technologies

CSols RFTrackIT – Capabilities & Demo