



laboratory systems integration

## CSols lab integration software reduces costs, turnaround time and transcription errors at Chemtest



Figure 1: New modern laboratories



Figure 2: New modern laboratories

Chemtest is a leading supplier of high quality, fast response analytical services where customer focus and continuous investment in personnel, processes and technology are at the forefront of their business strategy. Chemtest's laboratories are housed in modern, purposebuilt and secure accommodation in Newmarket. Here 65 highly qualified scientists carry out in excess of 6 million results determinations per year using a wide range of contemporary analytical instruments. CSols Ltd. has provided a number of software products and implementation services for the laboratory. These have had a significant effect upon Chemtest's activities, improving automation and hence reducing costs, average turnaround times and transcription errors. Furthermore a reduction in paper transcription and easily accessible audits of result processing using CSols software have helped Chemtest maintain their ISO/IEC 17025 and MCERTS accreditation and to successfully pass a number of audits.

#### Analysis at Chemtest

Analysis is undertaken at Chemtest for many industry sectors including Environmental, Geotechnical, Rail and Construction, Airports, Local Government as well as for Asbestos Management. Chemtest offers both fast and standard turnarounds to meet with clients' requirements and deadlines. The provision of sample containers and custody documentation, sample



Figure 3: New modern laboratories

collection and reporting costs are included in the price of analysis. Sample collection





laboratory systems integration

is provided utilising couriers and in-house transport to ensure timely and safe sample delivery to the laboratory so that samples are received in optimum condition for analysis with traceability from the sample origin. Chemtest also provides bespoke reporting formats in several different media to meet their clients' specifications.

#### **CSols & Compliance**

CSols software is helping Chemtest's laboratory workflow from sample reception and processing all the way to results delivery. In particular the following applications are used :

- PrepLab ensures samples are correctly labelled
- Links for LIMS helps analysts easily synchronise and process instrument results
- AqcTools is providing rapid feedback on the performance of each instrument

All the CSols software fully compliments the PerkinElmer LABWORKS LIMS and Totalchrom chromatography data handling system to help provide a complete system that allows Chemtest to easily meet its regulatory obligations for compliance for audit trailed sample and results histories for a wide range of methods. The following is a selection of methods where instruments (over 60 at the last count) use the CSols Links for LIMS application:

Photometric tests	(Discrete analyser)	PAH	(GC-MS)
Cyanide	(Continuous flow)	Phenols	(HPLC)
Metals	(ICP-MS)	TPH	(GC-FID)
Metals	(ICP-OES)	SVOCs	(GC-MS)
Soils - Moisture - init	ial weights		
Soils - Moisture - dry	v weights	PCBs	(GC-MS/ECD)
Sulphur	(TotalChrom)	VOCs	(GC-MS)
TOCS	(SciLabTOC)	TOC	(Carbon analyser)

Chemtest also participates in the following external proficiency schemes and CSols **AqcTools** is assisting in processing, storage and monitoring of the wide range of AQC performance data.

- CONTEST
- AQUACHECK
- LEAP
- AIMS
- RTC MCERTS



ENVIRON	MENTAL MA	NAGEMENT
SYS	FEM CERTIF	FICATE
This is to Certify	that the Environmental M	Aanagement System of
	CHEMTEST LT. Depot Road Newmarket CB8 0AL United Kingdom	D
has been	assessed by Sira Certific and found to comply w	cation Service
	ISO 14001:2004	1
	for	
The provision of chemical c	analysis services for envi and leachates).	ronmental samples (soils, waters
Certificate No: Date of Initial Certification: Date of Issue/Reissue: Renewal Due:	104029 02 June 2010 02 June 2010 12 May 2013	linda C
This certificate is subject to the compo	ny maintaining its system to the requir	of standards, which will be monitored by Sita.
The use of this Certificate and the Sira O	ertification Mark are subject to the Reg	slations Applicable to Helders of Sira Cortificates.
Sir	a Certification Se	Convice
	Dartford, Kent, United Kingdom, DA	1.44
The celeficate encloses the property of it	to and that be returned when requested, it may of	to republications enough and without change.



#### Sample Registration

After delivery samples are registered into Chemtest's LIMS, a confirmation of receipt is sent to the client; container specific barcoded sample labels and a folder label are produced. Technical Administration verifies that the customer order, pricing, tests / test suites are correct. The individual sample containers and folders are then labelled.

# 

laboratory systems integration



Figure 4: Sample registration

#### Sample preparation

Sample preparation at Chemtest utilises CSols **PrepLab** software. **PrepLab** streamlines and automates sample preparation for the analytical instrumental analysis. **PrepLab** is particularly well suited to the handling of solid samples such as soils as well as ores, metals, powders, soil, coal etc.



Figure 5: Sample preparation



Figure 6: Sample preparation area

Like all CSols software, **PrepLab** can be configured. At Chemtest users make use of barcode readers and capture weights directly from balances they are using in the preparation process. Chemtest also use **PrepLab** to produce further labels which help route prepared sub-samples to the correct analytical laboratory.





laboratory systems integration



#### PrepLab - Simplicity of use

Analysts typically have limited time to prepare samples and then setup each instrument so **PrepLab** operation is straightforward and easily handles the different instrument methods that are used for different sample types. **PrepLab** also allows users to record the physical layout of samples in preparation trays that are used or for worklists / sequence queues for different instruments. Unlike LIMS, in **PrepLab** the worklist run order is not forced on the analyst and so adding urgent samples, blanks, standards and QCs to the instrument's pre-sequence file is easy and straightforward and samples are both added though keyboard entry and by scanning using a barcode reader.

Automatic communication with PerkinElmer Labworks LIMS is also configured so that LIMS 'knows' that the sample is in the **PrepLab** tray. Any QC or blanks that are added are consistently named for the current method so that any subsequent results collected are easily extracted and processed by **AqcTools**. Upon completion of a sample tray two copies of the worklist are then produced by **PrepLab** – a physical and electronic copy. The first is printed and kept with the tray of samples to guide any further sample workup and the physical placement in the appropriate instrument autosampler. The second copy automatically updates LIMS to create an instrument worklist/sequence file that is used in the next phase of analysis.

#### Instrument setup



Instrument integration with the CSols **Links for LIMS** software is extensively used at Chemtest. For over 60 instruments of different types Links for LIMS provides a powerful user friendly interface for analysts to carry out their routine work.

Setting up an instrument is made easy as LIMS already has an instrument setup file. Analyst's simply use the Links for LIMS Setup program to select a worklist file that is matched to the sample tray. **Links for LIMS** then displays the samples and QCs in a layout suited for that specific instrument. If this matches the paper worklist attached to the tray, the run is accepted and an instrument compatible run file is created containing all the sample details and if required further biographical data, weights and dilution factors as required by the Standard Operating Procedure (SOP).





laboratory systems integration

🛿 Setup R	un - 02506	4313 o	n Kone								C
ile <u>E</u> dit <u>O</u>	utput <u>C</u> onfig	uration	Help								
III   <b>3</b> *	× 0 🖻	2 3	19	l							
	Sample	AS	Dil	TON	Nitrite	Phosphate	Ammonia	Chloride	Silicate	Alkalinity ·	•
1	Check Std	1/1	1								
2	0372358	1/2	1								
3	AQC 10	1/3	1								
4	0372359	1/4	1								
5	0372360	1/5	1								
6	0372362	1/6	1								
7	0372363	1/7	1								
8	0372365	1/8	1								
9	0372367	1/9	1								
10	0372368	1/ 10	1								
11	0372370	1/ 11	1								
12	0372371	1/ 12	1								
13	0372372	1/ 13	1								
14	0372373	1/ 14	1								
15	0372374	2/1	1								
16	0372378	2/2	1								
17	0372379	2/3	1								
18	0372381	2/4	1								
19	0372382	2/5	1								
20	0372383	2/6	1								
21	0372384	2/7	1								
22	0372385	2/8	1								
23	0372386	2/9	1								
24	0372389	2/10	1								
25	0372390	2/11	1								
26	0372392	2/ 12	1								
27	0372393	2/13	1								
28	0372394	2/14	1								
29	0372395	3/1	1								
30	AQC	3/2	1								
31	0372396	3/3	1								
32	0372397	3/4	1			-					
33	0372398	3/5	1								•
System Ma	essage>										
Samples	: 84	Method	: lo <del>w</del>	Access :	EDIT		V	iew : Works	heet		

Figure 7: Setting up a run for a Thermo/Kone Aquakem analyser showing the run list, autosampler position, dilution factor and range of results expected.

Depending upon the instrument system used, this run file can be opened up directly by the instrument software and used to run the instrument without further transcription.

#### Instrument results processing

Once an instrument has run, Chemtest's analysts use the **Links for LIMS** Result Review program to check their results before sending them to LIMS for final reporting. This gives them the chance to approve all of the data. The Result Reporting program is configured then to not only read the instrument data file but also to apply Chemtest's SOP settings for result treatment. This includes rules for rounding, treatment of blanks and AQCs, as well any calculations and handling of duplicates. AQC rules are applied with traffic light colour coding to highlight any problems which can be addressed easily by re-assigning samples for reanalysis in the next run or by discarding results (zapping) so that they are not sent to LIMS. Finally all acceptable results are transferred to LIMS by clicking a single function key in the software.





laboratory systems integration

Sample Name	Weight	Volume	Ag (Raw)	Ag	AI	B	Ba	Ca	Cd	Cr	Cu	Fe	К	Mg	Mn	Mo	Na	Ni	P
STD 1	1				20000						20000		20000					20000	
STD 2	1	XXXXXX	X0000X	XXXXXXX	X0000X	XXXXXXX	XXXXXXX	XXXXXXX	XXXXXXX	XXXXXXX	XXXXXXX								
STD 3	1	XXXXXX	X0000X	X0000X	XXXXXXX	200000	X0000X	X0000X	X0000X	X0000X	X0000X	X0000X	XXXXXXX	X0000X	XXXXXXX	XXXXXX	XXXXXXX	XXXXXXX	20
STD 4	1	XXXXXXX	xxxxxx	X0000X	X0000X	XXXXXXX	X0000X	XXXXXX	XXXXXX	X0000X	X0000X	x							
MULTINQC1	1		50	50.00	0.30	1.30	0.04	0.05	0.04	0.05	0.34	0.05	-0.10	0.04	0.05	0.08	0.21	0.04	
MULTINQC2	1		0.23	0.23	50.09	49.99	50.04	49.50	50.56	49.78	49.54	49.77	50.02	49.97	49.93	49.69	49.96	50.17	4
SAMPLE002	1	1.0	48	48.37	0.30	1.09	0.04	0.03	0.05	0.05	0.28	0.05	0.08	0.04	0.04	0.08	0.17	0.06	
SAMPLE004	1	1.0	0.16	0.16	50.03	49.85	49.94	49.18	50.54	49.47	49.34	49.43	49.89	49.75	49.64	50.51	49.79	50.17	4
WASH BLANK 1	1	10			0.19	8.56			0.19	0.23				0.16		0.59	0.61		
BLANK (X10)	1	10	0.48	0.48	0.27	5.16				0.10	1.77	0.09	0.64	0.01			0.36		
ADS 67674 (X10)	1	10	0.091	0.09	0.11	3.96	0.05	349.11	0.15	0.10	1.43	0.07	2.45	7.59	0.05	<0.00	12.40	0.17	
ADS 67675 (X10)	1	10	0.062	0.06	0.34	3.23	0.15	351.03	0.11	0.06	1.14	1.55	2.86	11.41	0.05	<0.00	14.08	0.10	<
ADS 67676 (X10)	1	10	0.031	0.03	0.33	2.62	0.03	399.33	0.11	0.07	0.91	0.09	4.26	17.78	0.27	<0.00	28.43	0.04	
WASH BLANK 2	1	10	0.94	0.94	0.46	2.13	0.06	0.30	0.08	0.08	0.76	0.06	<0.00	0.01	0.05	0.03	0.16	0.10	
MULTINGC1	1		0.085	0.00	50.21	49.60	0.00	40.20	50.77	49.61	49.37	10.01	50.14	40.04	49.74	50.55	49.97	50.21	
WASH BLANK 3	1	10	0.003	0.79	0.30	6.41	0.18	0.24	0.26	0.22	2.56	0.25	0.41	0.14	0.19	0.70	0.44	0.28	- 1
ADS 67677 (X10)	1	10	-0.018	<0.00	0.36	3.88	0.04	387.26	0.10	0.14	1.42	0.94	3.19	3.10	0.06	<0.00	11.41	0.16	-
ADS 67678 (X10)	1	10	0.027	0.03	0.26	2.74	0.04	274.80	0.07	0.06	0.98	0.05	1.34	0.01	0.03	<0.00	0.30	0.11	
ADS 67679 (X10)	1	10	-0.009	<0.00	0.20	2.22	<0.00	278.86	0.09	0.10	0.82	0.03	1.13	<0.00	0.03	<0.00	0.30	0.08	
ADS 67680 (X10)	1	10	0.076	0.08	0.05	1.81	0.19	338.46	<0.00	0.04	0.61	0.13	1.09	0.12	0.02	<0.00	20.19	0.04	
ADS 67681 (X10)	1	10	0.013	0.01	0.10	1.54	0.19	335.77	0.03	0.03	0.51	0.12	0.81	0.16	0.03	<0.00	28.61	0.05	
ADS 67682 (X10)	1	10	0.048	0.05	0.10	1.34	0.37	367.08	<0.00	0.06	0.47	0.16	0.99	0.14	0.03	<0.00	20.44	0.08	
WASH BLANK 4	1	10		0.87			0.06		0.02			0.06		0.02	0.04		<0.00		
MULTINQC1	1		49	49.00	0.24	0.29	0.00	0.01	0.01	0.01	0.03	0.01	-0.04	0.00	0.00	0.00	0.11	0.01	
MULTINQC2	1		0.063	0.06	50.23	49.83	50.07	49.89	51.20	49.97	49.87	49.77	50.45	50.54	50.02	50.85	50.09	50.22	5
WASH BLANK 5	1	10	0.74	0.74	0.32	5.94	0.25	0.32	0.23	0.28	2.58	0.28	<0.00	0.20	0.24	0.79	0.36	0.30	
ADS 67682B (X10)	1	10	0.049	0.05	0.14	3.24	0.35	362.96	0.08	0.12	1.19	0.21	1.79	0.15	0.06	<0.00	20.59	0.11	<
ADS 67683 (X10)	1	10	0.062	0.06	0.22	2.32	0.45	369.39	0.11	0.11	0.88	0.16	1.45	0.13	0.04	<0.00	18.68	0.08	
ADS 67684 (X10)	1	10	0.035	0.03	0.19	1.78	0.30	261.68	0.08	0.06	0.67	0.08	0.69	0.07	0.03	<0.00	19.79	0.08	
ADS 6/685 (X10)	1	10	-0.016	<0.00	0.02	1.50	0.30	266.50	0.03	0.06	0.58	0.08	<0.00	0.10	0.03	0.02	15.76	0.00	4
AUS 67665B (X10)	6	10	-0.004	×0.00	0.01	0.21	0.20	205.07	0.07	0.05	0.52	0.04	1.22	0.00	0.02	<0.00	15.79	0.01	
WASH BLANK 6	5	10	0.17	0.04	0.03	0.21	0.00	0.03	<0.00	0.00	0.04	0.00	<0.00	0.00	0.00	0.01	<0.00	0.01	
MULTINOC1	1	1.0	50	50.00	0.27	0.28	0.01	0.03	0.01	0.01	0.03	0.01	-0.10	0.00	0.00	0.00	0.10	0.00	
MULTINQC2	1		0.044	0.04	50.02	49.48	50.40	49.94	51.19	49.78	49.61	50.01	49.70	49.94	49.84	50.29	49.42	50,46	5
WASH BLANK 7	5	10	0.17	0.17	0.06	1.05	0.02	0.03	0.02	0.02	0.46	0.03	< 0.00	0.02	0.02	0.10	0.04	0.01	<
ADS 67674 (X2)	5	10	0.006	0.01	0.03	0.61	0.00	356.90	0.01	0.01	0.24	0.02	0.25	7.75	0.01	<0.00	11.96	<0.00	<
ADS 67675 (X2)	5	10	0.007	0.01	0.30	0.46	0.07	359.26	0.00	0.00	0.18	1.51	0.37	11.63	0.01	<0.00	13.72	0.01	<

Figure 9: Results reported from a PerkinElmer Optima showing colour coding of AQC samples

Both the electronic instrument setup and results uploading save Chemtest a considerable amount of time compared to doing this by hand. In addition the consistent interpretation of results provided by the Links for LIMS stored method has also saved time in reducing the amount of retesting.



#### AqcTools for routine analysis of AQCs

Chemtest is also currently rolling out another CSols software package called **AqcTools.** This will further help Chemtest in the

following ways :

 Immediate charting of the instrument AQC results within the Links for LIMS program to help spot problems with performance with an automatic upload to a central AqcTools database upon LIMS approval.





laboratory systems integration

- Online calculations of those specific variables and key parameters required by the accreditation bodies , for example, rolling bias, precision and so on.
- Central Database storage of the entire AQC dataset for easy processing and reporting.
- Ability to interact with discrete AQC values and flag results with comments and interpretation for later internal or external auditing.



Figure 10: Chart of results for each instrument run showing upper & lower control and action limits for the results dataset.

**AqcTools** is designed to run both within Links for LIMS as well as any laboratory workstation. This makes sure that Chemtest analysts get a consistent view of the data at all times. As there is no need to manually transfer results there are time and quality savings. As all results are stored in a full database environment, reporting is simplified saving a considerable amount of time, especially when supporting audits.





laboratory systems integration

#### Summary

CSols have over 21 years experience in producing solutions that assist manufacturing, services and process control labs, interfacing with a wide range of instrument types. CSols software tools have helped maximize Chemtest's investment across their entire laboratory workflow. In particular CSols solutions benefit the laboratory in the following ways :

- an easy to follow sample preparation system for capturing sample information and for labelling samples using PrepLab, which cuts out a large amount of manual transcription and the opportunity for identification error
- easier instrument processing by the 'earlier' creation of instrument run files including useful parameters such as dilution factor and sample weights saving valuable analyst time
- 'electronically' applies Chemtest SOPs to every instrument result to ensure the consistent interpretation of each result as well as the process the result must follow after review
- more responsive AQC system giving analysts an opportunity to identify problems as they occur to reduce wasteful analysis and instrument down time
- rapid data entry from instrument and central storage of all results in a database environment for easier reporting.

**Malcolm Avis, Technical Director at Chemtest** said "I have always found CSols very responsive to our needs. Our productivity gains made possible by their software, particularly with respect to sample preparation and instrument interfacing, have allowed us to reassign valuable resources to increase our throughput and reduce our turnaround times."

**Dr. Phil Goddard, CSols CEO** said "We constantly strive to achieve a close working relationship working with our customers, taking a hands-on approach, to ensure optimum efficiencies and performance. We believe that these are essential factors for creating and maintaining our customers' competitive advantage". We have enjoyed working with Chemtest and are delighted we have contributed to their continued success"





laboratory systems integration

#### **Further Information**

**Chemtest Limited**, Depot Road, Newmarket, CB8 0AL Tel: 01638 606070 Fax: 01638 606071 Web: <u>www.chemtest.co.uk</u> email : <u>malcolm.avis@chemtest.co.uk</u> Registered in England and Wales. Registration Number: 6511736 VAT Number: 927 131830

**CSols Ltd**, The Heath, Runcorn, Cheshire, WA7 4QX, UK Tel: +44 (0) 1928 513535 Fax: +44 (0) 7006 061106 Web: www.csols.com email: <u>chemtest@csols.com</u> Registered in England and Wales. Registration Number 3918227. VAT Number 539 693395

